

NEUROSENSE360

Complete Features & User Flow Documentation

EEG-Based Brain Wellness SaaS Platform

Version 2.0 | March 2026

Limitless Brain Lab (LBW)
<https://www.neurosense360.site>

Table of Contents

1. Platform Overview
2. Technology Stack
3. All Routes & Pages
4. Super Admin Features
5. Clinic Admin Features
6. Patient Features
7. Subscription Tiers & Pricing
8. User Flows
9. QEEG Report Generation Pipeline
10. Payment & Subscription Flow
11. Integrations & External Services
12. Database Schema
13. Landing Page Sections
14. Security & Compliance

1. Platform Overview

NeuroSense360 is a comprehensive, multi-role SaaS platform for EEG-based brain wellness assessment, clinic management, and patient engagement. It combines a public-facing marketing site with role-based dashboards for **Super Admins**, **Clinic Administrators**, and **Patients**.

The platform enables neurofeedback clinics to upload raw QEEG (quantitative EEG) data, process it through AI-powered algorithms, and generate branded brain wellness reports—all with built-in subscription management and payment processing.

Key Metrics Displayed

- 25,000+ members across the platform
- 2,500+ registered clinics
- Multi-currency payment support (8 currencies)
- AI-powered report generation in ~8-10 minutes

2. Technology Stack

LAYER	TECHNOLOGY
Frontend	React 18.2 + Vite 7.3 + Tailwind CSS 3.3
Routing	React Router v6
State / Forms	React Hook Form, React Query (TanStack)
UI	Lucide Icons, Framer Motion, Recharts, React Hot Toast
Backend	Node.js + Express.js (port 5000)
Database	Supabase (PostgreSQL + Auth + Storage)
AI Services	Google Gemini AI (primary), OpenAI (backup)
EEG Processing	qEEG Pro API + NeuroSense Cloud algorithms
Payments	Stripe (international) + Razorpay (India)
Email	Nodemailer via Gmail SMTP
PDF Generation	PDFKit + pdf-lib (13 modular page generators)
File Upload	Multer
Deployment	Vercel (frontend) + Render (backend)

3. All Routes & Pages

PUBLIC

Public Routes (No Authentication)

ROUTE	PURPOSE
/	Main landing page (18 sections, marketing site)
/about-us	Company information and team profiles
/pricing	Subscription tier comparison and pricing
/science	Scientific research and methodology
/guide-to-brainwaves	Educational content on brain waves
/technicians	Technician information and resources
/locations	Clinic locator with map integration
/contact	Contact form (sends to sweta.adata@gmail.com)
/neurosense-booking	Appointment booking system
/coach-certification	Coach certification program info
/100x-program-join	100x program enrollment form
/professional-onboarding	Professional registration form
/login	User login (all roles)
/register	User registration (patient or clinic)
/forgot-password	Password recovery
/reset-password	Password reset via token
/activation-pending	Account activation waiting page
/page/:slug	Dynamic CMS-managed static pages
/privacy-policy	Privacy policy

Wellness Protocol Pages (Public)

ROUTE	PURPOSE
-------	---------

/movers-section	MOVERS movement-based protocol
/ans-reset-protocol	Autonomic nervous system reset
/frequencies-music	Therapeutic frequencies and music
/supplements-nootropics	Supplement recommendations
/five-pillars	Five Pillars wellness framework
/brain-coach	Brain coaching program
/home-neurofeedback	Home neurofeedback protocol
/interactive-brain	Interactive brain visualization
/wallet	Patient wallet / credits

LBW (Limitless Brain Wellness) Pages

ROUTE	PURPOSE
/lbw	LBW main landing page
/lbw-updates	Project updates
/assessments	Brain wellness assessments
/coaching	Coaching programs
/lbw/dashboard	LBW user dashboard
/lbw/progress	Progress tracking
/lbw/content	Educational content library

SUPER ADMIN Admin Routes

ROUTE	FEATURE
/admin	System-wide dashboard and analytics
/admin/clinics	Clinic management (CRUD, approval, trials)
/admin/reports	All patient reports viewer
/admin/payments	Revenue tracking and transactions
/admin/analytics	System usage statistics
/admin/advanced-analytics	Cohort analysis, retention, churn
/admin/algorithm-processor	QEEG algorithm job processing
/admin/coaches	Coach management and assignments
/admin/alerts	Alert configuration and history
/admin/assessments	Assessment template management
/admin/static-pages	CMS - static page editor
/admin/settings	System configuration

CLINIC Clinic Routes

ROUTE	FEATURE
/clinic	Clinic overview dashboard
/clinic/patients	Patient management (CRUD, search, import)
/clinic/reports	Clinic's report viewer and downloads
/clinic/subscription	Purchase report packages
/clinic/usage	Usage analytics and trends
/clinic/admin	Clinic admin panel
/clinic/settings	Clinic configuration

PATIENT Patient Routes

ROUTE	FEATURE
-------	---------

/patient	Patient dashboard
/patient-dashboard	Patient dashboard (alternate)
/dashboard	Smart router (redirects by role)
/dashboard/*	Patient SPA sub-routes
/dashboard/subscription	Tier upgrade and management

4. Super Admin Features

1 Dashboard

System-wide analytics: active clinics, total patients, total reports, monthly revenue, system health, quick navigation to all sections.

2 Clinic Management

Create/edit/deactivate clinics, manage trial periods, approve new registrations, monitor usage and payment status, onboarding workflow.

3 Patient Reports

View all reports across all clinics, filter by clinic/date/status, download PDFs, manual upload, report approval workflow.

4 Payment History

Track all transactions, revenue breakdown by clinic, subscription tier monitoring, invoice management, payment status updates.

5 Analytics Dashboard

Real-time metrics, user growth trends, revenue analytics, clinic usage statistics, feature adoption, custom report generation.

6 Advanced Analytics

Deep-dive cohort analysis, retention metrics, churn analysis, custom dashboards, data export functionality.

7 Algorithm Processor

Process QEEG data through algorithms, manage configurations, view processing logs, validate results, bulk processing.

8 Coach Management

Create/manage coaches, assign to clinics, manage schedules, track performance, certification and specialization tracking.

9 Assessment Management

Create custom assessments, manage templates, track usage, view results, assessment analytics.

10 Alert System

Configure usage thresholds, email notification management, alert templates, alert history and logs.

11 Static Page CMS

Edit website pages, add/remove pages, manage content, SEO metadata, page versioning.

12 System Settings

Global configuration, email settings, payment gateway config, API key management, rate limiting.

5. Clinic Admin Features

1 Clinic Dashboard

Clinic-specific analytics, patient count, reports used vs. allowed, subscription status, recent activities, quick actions.

2 Patient Management

Add/edit patients, bulk import, search & filter, medical history, improvement focus areas, contact details, external ID tracking.

3 Report Viewer

View clinic's reports, download as PDF/EDF, filter by patient/date, report analytics per patient, progress timeline view.

4 Subscription Manager

View current tier, purchase additional report packages (per-10 model), billing history, upgrade/downgrade, trial tracking.

5 Usage Analytics

Reports used vs. allowed, monthly trends, patient activity, report processing times, storage usage, cost per report.

6 Clinic Settings

Edit clinic info, manage contact details, upload logo, configure address, team members, API credentials.

6. Patient Features

1 Patient Dashboard

Personal profile, brain assessment results, progress tracking, appointment booking, report downloads, recommended actions.

2 Brain Assessments

Brain Fitness Score, Brain Burnout Score, Brain Performance Score, Neuro Age Estimator, CogniFit Assessment, NeuroSense QEEG.

3 Progress Tracking

Historical data comparison, trend visualization, milestone achievements, weekly progress reports, goal tracking.

4 Wellness Programs

ANS Reset, MOVERS, Five Pillars, Frequencies & Music, Meditations, Supplements, Brain Coach, Home Neurofeedback (tier-gated).

5 Bookings

Schedule appointments, view available time slots, booking confirmations, reminders, cancellation/rescheduling.

6 Subscription

View/upgrade tier, compare plans, manage payment methods, billing history, instant feature unlock on upgrade.

7. Subscription Tiers & Pricing

FEATURE	FREE	BASIC ₹999	PRO ₹1,999	PREMIUM ₹3,499	ENTERPRISE ₹5,999
Profile & Brain Info	✓	✓	✓	✓	✓
Daily Inspiration	✓	✓	✓	✓	✓
ANS Reset Protocol	—	✓	✓	✓	✓
MOVERS Protocol	—	✓	✓	✓	✓
Five Pillars Framework	—	✓	✓	✓	✓
Frequencies & Music	—	—	✓	✓	✓
Meditations	—	—	✓	✓	✓
Supplements & Nootropics	—	—	✓	✓	✓
Brain Coach (1-on-1)	—	—	—	✓	✓
Home Neurofeedback	—	—	—	✓	✓
Assessments	—	—	—	✓	✓
Reports Included	0	10	25	50	100+

Multi-Currency Pricing

TIER	USD	INR	EUR	GBP	AED
Basic	\$9.99	₹849	€9.20	£7.90	AED 36.70
Pro	\$19.99	₹1,699	€18.39	£15.79	AED 73.37
Premium	\$34.99	₹2,999	€32.19	£27.59	AED 128.31
Enterprise	Custom pricing				

8. User Flows

Flow 1: Public Visitor → Patient

```
Landing Page → Click "Get Started"
  → Register (name, email, password, phone, address)
  → Email confirmation
  → Login → Patient Dashboard
  → Complete profile onboarding
  → Take initial brain assessments
  → View brain health scores
  → Access FREE tier content
  → Upgrade to BASIC/PRO/PREMIUM for more features
```

Flow 2: Public Visitor → Clinic Registration

```
Landing Page → Click "Clinic Login" / "Get Started"
  → Register as clinic (name, contact person, email, phone)
  → 30-day trial starts (50 free reports)
  → Clinic Dashboard → Onboarding wizard
  → Add first patient
  → Upload QEEG PDFs (Eyes Open + Eyes Closed)
  → AI pipeline generates report (~8-10 min)
  → Patient gets notified, views report
  → Monitor usage → Purchase more reports when low
```

Flow 3: Patient Report Viewing

```
Patient Login → Dashboard
  → Notification: "New report available"
  → View generated PDF report
  → See brain wave analysis & scores
  → Review personalized care plan
  → Track improvement metrics over time
  → Schedule follow-up appointment
```

Flow 4: Super Admin Daily Operations

Admin Login (superadmin@neuro360.com)

- Dashboard: check system analytics
- Review pending clinic registrations → Approve/reject
- Monitor revenue in Payment History
- Check report quality in Reports section
- Review algorithm processing queue
- Manage coach assignments
- Update CMS content if needed
- Review alerts for system issues

Flow 5: Payment / Subscription Upgrade

User hits usage limit or wants more features

- Navigate to Subscription page
- View current usage & available packages
- Select desired package/tier

IF location = India:

- Razorpay checkout (INR, UPI, cards, wallets)

ELSE:

- Stripe checkout (USD/EUR/GBP/AED + cards)
- Payment processed
- Instant quota increase / feature unlock
- Confirmation email + invoice generated

9. QEEG Report Generation Pipeline

1 Upload Phase

Clinic admin uploads Eyes Open PDF + Eyes Closed PDF. Files stored in temporary server location via Multer.

2 Logo Modification

Original EEG vendor logo removed from PDFs. NeuroSense branding applied. Modified PDFs re-saved with metadata preserved.

3 Cloud Upload

Both modified PDFs uploaded to Supabase qeeg-uploads storage bucket. Public access URLs generated for downstream processing.

4 qEEG Pro Processing

PDFs parsed to extract EEG data. Brain frequency bands calculated. Brain asymmetry analyzed. Coherence metrics computed. Standardized qEEG report generated.

5 NeuroSense Algorithm Analysis

Proprietary algorithms analyze brain wave patterns. Neurological risk assessment performed. Brain health scores calculated. Standardized report generated. Personalized care plan created.

6 Gemini AI Report Generation

Clinical insights synthesized from all data. Narrative report generated with recommendations. Care protocols personalized. Content formatted for PDF. NeuroSense branding applied.

7 Report Finalization

Final PDF stored in Supabase neurosense-reports bucket. Patient record updated with report link. Audit log entry created. Patient notification sent via email. Reports-used counter incremented for clinic.

Estimated Total Processing Time: 8–10 minutes from upload to final report delivery.

PDF Report Structure (13 Modules)

#	MODULE	CONTENT
1	coverPage.js	Branded cover with patient info and date
2	introductionPage.js	Introduction to brain wellness assessment
3	brainTypePage.js	Brain type classification and description
4	brainwaveProfilesPage.js	Delta, Theta, Alpha, Beta, Gamma analysis
5	brainTopology.js	Brain region mapping and topology
6	brainMapComparisonPage.js	Before/after or norm comparison maps
7	cognitionPage.js	Cognitive function assessment
8	yourNumbersPage.js	Detailed numerical scores and metrics
9	numbersAtGlance.js	Summary scorecard view
10	metricsPages.js	Detailed metrics with graphs
11	congratulationsPage.js	Summary and next steps
12	appendix.js	Technical appendix and references
13	staticPages.js	Legal disclaimers and static info

10. Payment & Subscription Flow

Clinic Subscription Model

- 1. **Registration:** 30-day trial period with 50 free reports, no payment required
- 2. **Trial Warnings:** Email alerts at 7 days, 3 days, and 1 day before expiration
- 3. **Purchase:** Choose package (10 / 25 / 50 / 100 reports) with volume discounts
- 4. **Payment:** India → Razorpay (INR, UPI, cards) | International → Stripe (multi-currency)
- 5. **Activation:** Instant quota increase, confirmation email, invoice generated

Patient Subscription Model

- 1. **Free Tier:** Basic brain info, profile, limited assessments, 0 reports
- 2. **Upgrade:** Click "Upgrade" in dashboard, see tier comparison
- 3. **Payment:** Currency auto-detected by location, Stripe or Razorpay
- 4. **Activation:** Features unlock instantly, confirmation email sent

Supported Payment Methods

GATEWAY	REGION	METHODS	CURRENCIES
Razorpay	India	Credit/Debit Cards, UPI, Net Banking, Wallets	INR
Stripe	International	Credit/Debit Cards, Apple Pay, Google Pay	USD, EUR, GBP, AED, AUD, CAD, SGD

11. Integrations & External Services

SERVICE	PURPOSE	USAGE
Supabase	Database, Auth, Storage	PostgreSQL for all data, Auth for user management, 3 storage buckets (neurosense-reports, clinic-logos, qeeg-uploads)
Google Gemini AI	Report Generation	PDF report content, care plan synthesis, personalized recommendations. Models: gemini-3-pro, flash, pro
OpenAI	Backup AI	Fallback for report generation when Gemini is unavailable

Stripe	International Payments	Subscription checkout, one-time purchases, 8 currencies, invoices
Razorpay	Indian Payments	INR transactions, UPI, subscription management, payment links
qEEG Pro	EEG Processing	EDF file processing, brain wave analysis, quantitative EEG computation
NeuroSense Cloud	Algorithm Processing	Proprietary brain wave algorithms, risk assessment, care plan generation
Gmail SMTP	Email	Transactional emails, notifications, payment confirmations, appointment reminders
Google Maps	Location	Clinic location detection, address geocoding, map display

12. Database Schema

TABLE	PURPOSE	KEY FIELDS
clinics	Clinic information	name, email, phone, address, is_active, trial_end_date, reports_allowed
patients	Patient records	name, email, phone, dob, gender, clinic_id, external_id, medical_history
profiles	User profiles (auth)	user_id, role, clinic_id, display_name
reports	Report metadata	patient_id, clinic_id, pdf_url, status, created_at, type
appointments	Booking data	patient_id, clinic_id, date, time_slot, status, type
subscriptions	Subscription tracking	user_id, tier, status, start_date, end_date, payment_method
payment_history	Transaction records	user_id, amount, currency, gateway, status, invoice_id
coaches	Coach information	name, specialization, clinic_ids, schedule, certification
algorithm_results	QEEG processing results	patient_id, report_id, scores, brain_type, metrics
access_logs	Feature access tracking	user_id, feature, timestamp, ip_address
clinical_forms	Assessment form data	patient_id, form_type, responses, score, completed_at

Supabase Storage Buckets

BUCKET	PURPOSE	ACCESS
neurosense-reports	Final generated PDF reports for patients	Authenticated (RLS)
clinic-logos	Clinic branding and logo images	Public
qeeg-uploads	Raw QEEG/EDF files uploaded by clinics	Authenticated (RLS)

13. Landing Page Sections (18 Total)

#	SECTION	DESCRIPTION
1	Announcement Bar	Limited-time promotional offers
2	Scroll Progress	Visual progress indicator while scrolling
3	Navigation Bar	Fixed top menu with CTAs
4	Hero Section	Bold headline, social proof stats, primary CTAs
5	How It Works	4-step process visualization
6	Stats Section	7 animated statistics counters
7	Benefits	6 use-case cards (stress, sleep, focus, ADHD, etc.)
8	Press & Media	6 publication logos and mentions
9	Testimonials	3 customer stories with photos
10	Trust Badges	Partner hospitals and certifications
11	What's Included	Feature breakdown by tier
12	Comparison Table	NeuroSense vs. Medication vs. Therapy
13	Expert Team	Team member profiles and credentials
14	Results Timeline	5-phase improvement journey
15	Pricing	Subscription tier cards
16	FAQ Accordion	8 commonly asked questions
17	Final CTA	Strong closing call-to-action
18	Footer	Links, newsletter signup, social media

14. Security & Compliance

AREA	IMPLEMENTATION
------	----------------

Authentication	Password hashing with bcryptjs, JWT tokens, session cookies, Supabase Auth
Authorization	Role-based access control (RBAC): super_admin, clinic_admin, patient. ProtectedRoute component enforces access.
Data Protection	Supabase Row Level Security (RLS), HTTPS encryption, secure password recovery flow
Audit Trail	Access logs for sensitive data, report download tracking, payment transaction logging, admin action logging
Email Security	Email verification before account activation, Gmail App Passwords (not plain passwords)
Compliance	HIPAA considerations (medical data), GDPR data handling, privacy policy page, terms of service