

INTRODUCTION


The qEEG report provided by NeuroSense EEG is intended for informational, educational, and wellness purposes only. It is designed to help individuals and neurofeedback professionals better understand brainwave patterns and to support decisions related to neurofeedback training for non-medical cognitive enhancement. This report is not intended to diagnose, treat, cure, mitigate, or prevent any medical condition, and it should not be used as a substitute for consultation with a licensed healthcare provider.

Users are encouraged to consult with a qualified healthcare provider regarding any medical concerns or before starting any new treatment or therapy. NeuroSense EEG's qEEG analysis application is not a replacement for the individualized care provided by medical professionals.

EEG RECORDING

The 10-20 system is the internationally recognized method used for placing electrodes on the scalp during an EEG (Electroencephalogram) recording. It is named for the standardized distances between electrode positions, which are either 10% or 20% of the total front-to-back or right-to-left measurement of the head. This system ensures consistent and reproducible electrode placement, allowing for accurate brainwave measurement across individuals.

Electrodes are positioned over specific areas of the brain, corresponding to functional regions like the frontal, temporal, parietal, and occipital lobes, helping clinicians and researchers capture electrical activity associated with various cognitive and neurological functions. The 10-20 system is widely used in clinical diagnostics and research to assess brain activity related to conditions such as epilepsy, sleep disorders, and other neurological disorders.



A segment of raw EEG signal from the 19 electrode locations in the Eyes Closed condition.

1. Arns, M., de Ridder, S., Strehl, U., Breteler, M., & Coenen, A. (2009). Efficacy of neurofeedback treatment in ADHD: The effects on inattention, impulsivity, and hyperactivity: A meta-analysis. *Clinical EEG and Neuroscience*, 40(3), 180-189. doi: 10.1177/155005940904000311 <https://www.ncbi.nlm.nih.gov/pubmed/19715181>.
2. Micoulaud-Franchi, J.-A., Geoffroy, P. A., Fond, G., Lopez, R., Bioulac, S., Philip, P. (2014). EEG neurofeedback treatments in children with ADHD: An update meta-analysis of randomized controlled trials. *Frontiers in Human Neuroscience*, 8(906), 1-7. doi: 10.3389/fnhum.2014.00906 <https://www.ncbi.nlm.nih.gov/pubmed/25431555>.
3. Steiner, N. J., Frenette, E. C., Rene K. M., Brennan, R. T., & Perrin, E. C. (2014). In-school neurofeedback training for ADHD: Sustained improvements from a randomized control trial. *Pediatrics*, 133(3), 483-492. doi: 10.1542/peds.2013-2059. <https://www.ncbi.nlm.nih.gov/pubmed/24534402>.
4. Wigton, N. L., & Krigbaum, G. (2015). Attention, executive function, behavior, and electrocortical function, significantly improved with 19-channel z-score neurofeedback in a clinical setting: A pilot study. *Journal of Attention Disorders*, [e-pub ahead of print]. doi: 10.1177/1087054715577135, <https://www.ncbi.nlm.nih.gov/pubmed/25823743>.
5. Escolano, Carlos, and Javier Minguez. "EEG-Based Upper Alpha Neurofeedback Training Improves Working Memory Performance." *IEEE Elore*, 2011, <https://ieeexplore.ieee.org/abstract/document/6090651/>.
6. Gruzelier, John H. "EEG-Neurofeedback for Optimising Performance. III: A Review of Methodological and Theoretical Considerations." *Neuroscience & Biobehavioral Reviews*, Pergamon, 29 Mar. 2014, <https://www.sciencedirect.com/science/article/abs/pii/S0149763414000700>.
7. Woelfel, Benedikt, et al. "Neurofeedback Training of the Upper Alpha Frequency Band in EEG Improves Cognitive Performance." *NeuroImage*, Academic Press, 17 Sept. 2010, <https://www.sciencedirect.com/science/article/abs/pii/S105381191001181X>.